




AT135 rabbit pAb

Catalog No	YP-Ab-08733
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ATP13A5 UNQ488/PRO1004
Immunogen	Synthesized peptide derived from human AT135 AA range: 811-861
Specificity	This antibody detects endogenous levels of AT135 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AT135;ATP13A5 UNQ488/PRO1004
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	catalytic activity:ATP + H(2)O = ADP + phosphate.,sequence caution:Translation N-terminally extended.,sequence caution:Unlikely isoform. Cloning artifact., similarity:Belongs to the cation transport ATPase (P-type) family. Type V subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4°over night