




FA40B rabbit pAb

Catalog No	YP-Ab-08737
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM40B KIAA1170
Immunogen	Synthesized peptide derived from human FA40B AA range: 602-652
Specificity	This antibody detects endogenous levels of FA40B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA40B;FAM40B KIAA1170
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Enriched in lamellipodia.
Tissue Specificity	
Function	similarity:Belongs to the FAM40 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night