




EMC2 rabbit pAb

Catalog No	YP-Ab-08746
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	EMC2 KIAA0103 TTC35
Immunogen	Synthesized peptide derived from human EMC2 AA range: 103-153
Specificity	This antibody detects endogenous levels of EMC2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EMC2;ER membrane protein complex subunit 2;KIAA0103;TPR repeat protein 35;EMC2 KIAA0103 TTC35
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Peripheral membrane protein ; Cytoplasmic side . May also localize to the nuclear envelope. .
Tissue Specificity	
Function	similarity:Contains 3 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night