




C2G1P rabbit pAb

Catalog No	YP-Ab-08763
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CYP2G1P CYP2GP1
Immunogen	Synthesized peptide derived from human C2G1P AA range: 13-63
Specificity	This antibody detects endogenous levels of C2G1P at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2G1P;CYP2G1P CYP2GP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:Could be the product of a pseudogene. In contrast to other members of the family, it is shorter at the N-terminus.,cofactor:Heme group., similarity:Belongs to the cytochrome P450 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night