




NAA16 rabbit pAb

Catalog No	YP-Ab-08766
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NAA16 NARG1L NAT2
Immunogen	Synthesized peptide derived from human NAA16 AA range: 800-850
Specificity	This antibody detects endogenous levels of NAA16 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAA16;N-alpha-acetyltransferase 16;NatA auxiliary subunit;NARG1 like protein; NARG1-like protein;NAT2;NAA16 NARG1L NAT2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,transcription factor complex,cytoplasm,NatA complex,extracellular exosome,
Tissue Specificity	
Function	function:May belong to a complex displaying N-terminal acetyltransferase activity.,similarity:Contains 7 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night