




ZDH12 rabbit pAb

Catalog No	YP-Ab-08775
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ZDHHC12 ZNF400 PSEC0008
Immunogen	Synthesized peptide derived from human ZDH12 AA range: 83-133
Specificity	This antibody detects endogenous levels of ZDH12 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZDH12; ZDHHC12 ZNF400 PSEC0008
Observed Band	
Calculated Molecular Weight	29kD
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	catalytic activity: Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA., domain: The DHHC domain is required for palmitoyltransferase activity., similarity: Belongs to the DHHC palmitoyltransferase family., similarity: Contains 1 DHHC-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4°over night