




TBC12 rabbit pAb

Catalog No	YP-Ab-08783
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TBC1D12 KIAA0608
Immunogen	Synthesized peptide derived from human TBC12 AA range: 596-646
Specificity	This antibody detects endogenous levels of TBC12 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBC1 domain family member 12;TBC12_HUMAN;TBC1D12;TBC12;TBC1D12 KIAA0608
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endosome .
Tissue Specificity	
Function	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night