




CB057 rabbit pAb

Catalog No	YP-Ab-08785
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	C2orf57
Immunogen	Synthesized peptide derived from human CB057 AA range: 259-309
Specificity	This antibody detects endogenous levels of CB057 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CB057; C2orf57
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis. Detected in germ cells at all stages of the seminiferous epithelium, strong expression in elongating spermatids (at protein level) (PubMed:26168773).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night