




NUDT8 rabbit pAb

Catalog No	YP-Ab-08786
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NUDT8
Immunogen	Synthesized peptide derived from human NUDT8 AA range: 89-139
Specificity	This antibody detects endogenous levels of NUDT8 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.6.1.-;Mitochondrial coenzyme A diphosphatase NUDT8;Nucleoside diphosphate-linked moiety X motif 8;Nudix motif 8;NUDT8
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	cofactor:Magnesium or manganese.,function:Probably mediates the hydrolysis of some nucleoside diphosphate derivatives.,similarity:Belongs to the Nudix hydrolase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night