




YN3949

DIEXF rabbit pAb

Catalog No	YP-Ab-08791
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	DIEXF C1orf107 DEF
Immunogen	Synthesized peptide derived from human DIEXF AA range: 296-346
Specificity	This antibody detects endogenous levels of DIEXF at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf107;DEF;DIEXF;DIEXF C1orf107 DEF
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Expressed in colon.
Function	function:Regulates the p53 pathway to control the expansion growth of digestive organs.,similarity:Belongs to the def family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4°over night