




PRRP rabbit pAb

Catalog No	YP-Ab-08802
Isotype	IgG
Reactivity	Human;Rat
Applications	IHC;IF
Gene Name	PRLH PRH
Immunogen	Synthesized peptide derived from human PRRP AA range: 16-66
Specificity	This antibody detects endogenous levels of PRRP at Human/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PrRPR;PrRP receptor;hGR3;GR3;GPR10;PRRP;PRLH PRH
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Medulla oblongata and hypothalamus.
Function	function:Stimulates prolactin (PRL) release and regulates the expression of prolactin through its receptor GPR10. May stimulate lactotrophs directly to secrete PRL.,tissue specificity:Medulla oblongata and hypothalamus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).