




ARP5L rabbit pAb

Catalog No	YP-Ab-08803
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ARPC5L
Immunogen	Synthesized peptide derived from human ARP5L AA range: 4-54
Specificity	This antibody detects endogenous levels of ARP5L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARP5L; ARPC5L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	function: May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks., similarity: Belongs to the ARPC5 family., subunit: May be a component of the Arp2/3 complex in which it may replace ARPC5.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night