




P3IP1 rabbit pAb

Catalog No	YP-Ab-08805
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	PIK3IP1 HGFL
Immunogen	Synthesized peptide derived from human P3IP1 AA range: 107-157
Specificity	This antibody detects endogenous levels of P3IP1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HGFL; hHGFL(S); Kringle domain containing protein HGFL; Kringle domain-containing protein HGFL; MGC1733; P3IP1_HUMAN; Phosphoinositide 3 kinase interacting protein 1; Phosphoinositide-3-kinase-interacting protein 1; PIK3IP1; P3IP1; PIK3IP1 HGFL
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	similarity: Contains 1 kringle domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night