




YC002 rabbit pAb

Catalog No	YP-Ab-08814
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	UNQ6490/PRO21339
Immunogen	Synthesized peptide derived from human YC002 AA range: 66-116
Specificity	This antibody detects endogenous levels of YC002 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YC002;UNQ6490/PRO21339
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night