




TAD2B rabbit pAb

Catalog No	YP-Ab-08819
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TADA2B ADA2B
Immunogen	Synthesized peptide derived from human TAD2B AA range: 220-270
Specificity	This antibody detects endogenous levels of TAD2B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p68;GANDS;MRD18;P66beta;GATAD2B;TAD2B;TADA2B ADA2B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Coactivates PAX5-dependent transcription together with either SMARCA4 or GCN5L2.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 1 SANT domain.,similarity:Contains 1 ZZ-type zinc finger.,subunit:Interacts with GCN5L2, SMARCA4, SMARCE1 and PAX5. Component of the TFTC-HAT complex.,
Background	TADA2B functions as a transcriptional adaptor protein that potentiates transcription through coordination of histone acetyltransferase (HAT) activity and by linking activation factors to basal transcriptional machinery (Barlev et al., 2003 [PubMed 12972612]).[supplied by OMIM, Apr 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night