




K0284 rabbit pAb

Catalog No	YP-Ab-08820
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KIAA0284 FAM68C
Immunogen	Synthesized peptide derived from human K0284 AA range: 454-504
Specificity	This antibody detects endogenous levels of K0284 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K0284;KIAA0284 FAM68C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	function:Plays a role in microtubule organization.,similarity:Belongs to the CEP170 family.,similarity:Contains 1 FHA domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night