




DTWD1 rabbit pAb

Catalog No	YP-Ab-08828
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	DTWD1 MDS009
Immunogen	Synthesized peptide derived from human DTWD1 AA range: 31-81
Specificity	This antibody detects endogenous levels of DTWD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DTW domain containing 1;DTWD1;MDS009;DTWD1 MDS009
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	domain:Contains 1 DXTW motif.,similarity:Belongs to the DTW family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night