




RAB44 rabbit pAb

Catalog No	YP-Ab-08830
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RAB44
Immunogen	Synthesized peptide derived from human RAB44 AA range: 372-422
Specificity	This antibody detects endogenous levels of RAB44 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RASL13;RASD3;RAB44
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	
Function	similarity:Belongs to the small GTPase superfamily. Rab family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night