




KBTB4 rabbit pAb

Catalog No	YP-Ab-08833
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KBTBD4 BKLHD4
Immunogen	Synthesized peptide derived from human KBTB4 AA range: 393-443
Specificity	This antibody detects endogenous levels of KBTB4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KBTB4;KBTBD4 BKLHD4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,SCF ubiquitin ligase complex,
Tissue Specificity	
Function	similarity:Contains 1 BACK (BTB/Kelch associated) domain.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 Kelch repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night