




DCA15 rabbit pAb

Catalog No	YP-Ab-08841
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DCAF15 C19orf72
Immunogen	Synthesized peptide derived from human DCA15 AA range: 251-301
Specificity	This antibody detects endogenous levels of DCA15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DCA15;DCAF15 C19orf72
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function:May be involved in ubiquitination and degradation through a DBB1-CUL4 E3 protein-ubiquitin ligase.,pathway:Protein modification; protein ubiquitination.,subunit:Interacts with DDB1 and CUL4A.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night