



YN4012

ARFP2 rabbit pAb


Catalog No	YP-Ab-08852
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ARFIP2 POR1
Immunogen	Synthesized peptide derived from human ARFP2 AA range: 137-187
Specificity	This antibody detects endogenous levels of ARFP2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARFP2;ARFIP2 POR1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus . Golgi apparatus, trans-Golgi network membrane .
Tissue Specificity	
Function	function:Putative target protein of ADP-ribosylation factor. Involved in membrane ruffling.,similarity:Contains 1 AH domain.,subunit:Interacts with RAC1 by binding directly to it. Specifically binds to GTP-bound ARF1 and ARF6, but binds to RAC1.GTP and RAC1.GDP with similar affinities.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4°over night