




AN13B rabbit pAb

Catalog No	YP-Ab-08860
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ANKRD13B
Immunogen	Synthesized peptide derived from human AN13B AA range: 544-594
Specificity	This antibody detects endogenous levels of AN13B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AN13B;ANKRD13B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane . Late endosome . Early endosome . Interaction with EGFR may enhance association with the cell membrane.
Tissue Specificity	
Function	similarity:Contains 2 ANK repeats.,similarity:Contains 3 UIM (ubiquitin-interacting motif) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night