



YN4021

CLVS2 rabbit pAb


Catalog No	YP-Ab-08861
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	CLVS2 C6orf212 C6orf213 RLBP1L2
Immunogen	Synthesized peptide derived from human CLVS2 AA range: 26-76
Specificity	This antibody detects endogenous levels of CLVS2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C6orf212;C6orf213;clathrin vesicle-associated Sec14 protein 2;Clavesin 2; Clavesin-2;CLVS2;CLVS2 C6orf212 C6orf213 RLBP1L2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle . Early endosome membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	
Background	This gene encodes a protein that belongs to the SEC14/CRAL-TRIO family of proteins. A similar protein in rat is thought to function in the endosomal pathway between early endosomes and mature lysosomes. [provided by RefSeq, Aug 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night