



OCAD1 rabbit pAb


Catalog No	YP-Ab-08866
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	OCIAD1 OCIA
Immunogen	Synthesized peptide derived from human OCAD1 AA range: 114-164
Specificity	This antibody detects endogenous levels of OCAD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Asrij; FLJ2455; MGC11172; OCAD1_HUMAN; OCIA; OCIA domain containing 1; OCIA domain-containing protein 1; OCIAD1; Ovarian carcinoma immunoreactive antigen; TPA18; OCAD1; OCIAD1 OCIA
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endosome .
Tissue Specificity	Isoform 1 is highly expressed in many tissues, including testis, brain, placenta, ovary, prostate and mammary gland. Isoform 2 expression is restricted to the central nervous system including brain, cerebellum and spinal cord.
Function	similarity: Belongs to the OCIAD1 family., similarity: Contains 1 OCIA domain., tissue specificity: Isoform 1 is highly expressed in many tissues, including testis, brain, placenta, ovary, prostate and mammary gland. Isoform 2 expression is restricted to the central nervous system including brain, cerebellum and spinal cord.,.
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night