




UBFL1 rabbit pAb

Catalog No	YP-Ab-08869
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UBTFL1 C11orf27
Immunogen	Synthesized peptide derived from human UBFL1 AA range: 223-273
Specificity	This antibody detects endogenous levels of UBFL1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBFL1;UBTFL1 C11orf27
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Mainly cytoplasmic from the 4-cell stage to the morula stage. Becomes nuclear at the blastocyst stage. .
Tissue Specificity	
Function	
Background	UBTFL1 is a preimplantation-specific gene and is involved in early development, implantation, and embryonic stem (ES) cell derivation (summary by Yamada et al., 2010 [PubMed 19915186]).[supplied by OMIM, Jan 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night