




XKR8 rabbit pAb

Catalog No	YP-Ab-08880
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	XKR8 XRG8
Immunogen	Synthesized peptide derived from human XKR8 AA range: 111-161
Specificity	This antibody detects endogenous levels of XKR8 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	XKR8;XKR8 XRG8
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cytoplasm, perinuclear region .
Tissue Specificity	
Function	similarity:Belongs to the XK family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night