




MYADM rabbit pAb

Catalog No	YP-Ab-08883
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	MYADM UNQ553/PRO1110
Immunogen	Synthesized peptide derived from human MYADM AA range: 80-130
Specificity	This antibody detects endogenous levels of MYADM at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	myadm; MYADM_HUMAN; myeloid associated differentiation marker; Myeloid upregulated protein; Myeloid-associated differentiation marker; Protein SB135; SB135; MYADM UNQ553/PRO1110
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Not detected in thymus.
Function	similarity: Belongs to the MAL family., similarity: Contains 2 MARVEL domains., tissue specificity: Widely expressed. Not detected in thymus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night