




FA96A rabbit pAb

Catalog No	YP-Ab-08897
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM96A
Immunogen	Synthesized peptide derived from human FA96A AA range: 58-108
Specificity	This antibody detects endogenous levels of FA96A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA96A;FAM96A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Substantially enriched in macrophages.
Function	similarity:Belongs to the UPF0195 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night