




IV4F8 rabbit pAb

Catalog No	YP-Ab-08904
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	IGHV4OR15-8 VSIG6
Immunogen	Synthesized peptide derived from human IV4F8 AA range: 36-86
Specificity	This antibody detects endogenous levels of IV4F8 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IV4F8;IGHV4OR15-8 VSIG6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	external side of plasma membrane,immunoglobulin complex, circulating,blood microparticle,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night