




DUPD1 rabbit pAb

Catalog No	YP-Ab-08916
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	DUPD1 DUSP27
Immunogen	Synthesized peptide derived from human DUPD1 AA range: 80-130
Specificity	This antibody detects endogenous levels of DUPD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DUPD1; Dual specificity phosphatase 29; DUSP27; DUSP29; EC:3.1.3.16; EC:3.1.3.48; DUPD1 DUSP27
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	catalytic activity: A phosphoprotein + H ₂ O = a protein + phosphate., catalytic activity: Protein tyrosine phosphate + H ₂ O = protein tyrosine + phosphate., similarity: Belongs to the protein-tyrosine phosphatase family. Non-receptor class dual specificity subfamily., similarity: Contains 1 tyrosine-protein phosphatase domain., tissue specificity: Skeletal muscle, liver and adipose tissue.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night