




TTL13 rabbit pAb

Catalog No	YP-Ab-08918
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TTL13
Immunogen	Synthesized peptide derived from human TTL13 AA range: 423-473
Specificity	This antibody detects endogenous levels of TTL13 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTL13;TTLL13
Observed Band	
Calculated Molecular Weight	
Cell Pathway	microtubule,
Tissue Specificity	
Function	function:Polyglutamylase which preferentially modifies alpha-tubulin. Involved in the side-chain elongation step of the polyglutamylation reaction rather than in the initiation step.,similarity:Belongs to the tubulin--tyrosine ligase family., similarity:Contains 1 TTL domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night