




ABH15 rabbit pAb

Catalog No	YP-Ab-08925
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ABHD15 UNQ6510/PRO21435
Immunogen	Synthesized peptide derived from human ABH15 AA range: 101-151
Specificity	This antibody detects endogenous levels of ABH15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABH15_HUMAN;Abhd15;Abhydrolase domain containing protein 15; Abhydrolase domain-containing protein 15;Alpha/beta hydrolase domain-containing protein 15;hypothetical protein LOC116236;LOC116236;Protein ABHD15;RGD134598;ABH15;ABHD15 UNQ6510/PRO21435
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Belongs to the AB hydrolase superfamily. AB hydrolase 4 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night