




TM218 rabbit pAb

Catalog No	YP-Ab-08926
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM218
Immunogen	Synthesized peptide derived from human TM218 AA range: 32-82
Specificity	This antibody detects endogenous levels of TM218 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM218; TMEM218
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Cell projection, cilium . Localizes at the transition zone, a region between the basal body and the ciliary axoneme. .
Tissue Specificity	
Function	similarity: Belongs to the TMEM218 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night