



SRPX rabbit pAb


Catalog No	YP-Ab-08937
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SRPX ETX1
Immunogen	Synthesized peptide derived from human SRPX AA range: 234-284
Specificity	This antibody detects endogenous levels of SRPX at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Down-regulated by v-SRC;DRS;DRS protein;Drs-1;Drs-2;Epididymis secretory sperm binding protein Li 83p;ETX1;EXT1;HEL S 83p;mDRS-1;mDRS-2;SrpX;SRPX_HUMAN;SRPX1;Sushi repeat containing protein;X chromosome;X linked;Sushi repeat-containing protein SRPX;Sushi-repeat-containing protein;Sushi-repeat-containing protein SRPX;SRPX ETX1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell surface . Possibly surface of photoreceptor cell.
Tissue Specificity	Retina and heart; less in placenta, pancreas, lung, liver, skeletal muscle, kidney and brain.
Function	alternative products:Additional isoforms seem to exist,function:May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.,similarity:Contains 1 HYR domain.,similarity:Contains 3 Sushi (CCP/SCR) domains.,subcellular location:Possibly surface of photoreceptor cell.,tissue specificity:Retina and heart; less in placenta, pancreas, lung, liver, skeletal muscle, kidney and brain.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night