




FRPD3 rabbit pAb

Catalog No	YP-Ab-08941
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FRMPD3 KIAA1817
Immunogen	Synthesized peptide derived from human FRPD3 AA range: 1451-1501
Specificity	This antibody detects endogenous levels of FRPD3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FRPD3;FRMPD3 KIAA1817
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoskeleton,
Tissue Specificity	
Function	similarity:Contains 1 FERM domain.,similarity:Contains 1 PDZ (DHR) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night