




TFPT rabbit pAb

Catalog No	YP-Ab-08959
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TFPT INO80F
Immunogen	Synthesized peptide derived from human TFPT AA range: 47-97
Specificity	This antibody detects endogenous levels of TFPT at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TCF3 fusion partner; Protein FB1; INO80F; INO80 complex subunit F; FB1; TFPT; TFPT INO80F
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	disease: A chromosomal aberration involving TFPT is a cause of pre-B-cell acute lymphoblastic leukemia (B-ALL). Inversion inv(19)(p13;q13) with TCF3., subunit: Interacts with NOL3 (By similarity). Component of the chromatin-remodeling INO80 complex, at least composed of ACTL6A, ACTR5, ACTR8, RVBL1, RVBL2, INO80, INO80B, INO80C, INO80D and INO80E.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night