




RPR1B rabbit pAb

Catalog No	YP-Ab-08961
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RPRD1B C20orf77 CREPT
Immunogen	Synthesized peptide derived from human RPR1B AA range: 248-298
Specificity	This antibody detects endogenous levels of RPR1B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RPR1B;RPRD1B C20orf77 CREPT
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Preferentially expressed in a range of tumor tissues including colon, lung, liver, breast, prostate, stomach, uterine endometrium and cervical cancers with higher levels in tumors than in adjacent non-tumor tissue (at protein level).
Function	sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the UPF0400 (RTT103) family.,similarity:Contains 1 CID domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4°over night