




FA73B rabbit pAb

Catalog No	YP-Ab-08962
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM73B C9orf54 PSEC0112
Immunogen	Synthesized peptide derived from human FA73B AA range: 6-56
Specificity	This antibody detects endogenous levels of FA73B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA73B;FAM73B C9orf54 PSEC0112
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Mitochondrion outer membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the FAM73 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night