




CD006 rabbit pAb

Catalog No	YP-Ab-08972
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	C4orf6 AC1
Immunogen	Synthesized peptide derived from human CD006 AA range: 6-56
Specificity	This antibody detects endogenous levels of CD006 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD006;C4orf6 AC1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in neuroblastoma.
Function	tissue specificity:Expressed in neuroblastoma.,
Background	This gene is expressed in neuroblastoma; however, the function of this gene is not yet determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).