




COBL1 rabbit pAb

Catalog No	YP-Ab-08973
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	COBLL1 KIAA0977
Immunogen	Synthesized peptide derived from human COBL1 AA range: 725-775
Specificity	This antibody detects endogenous levels of COBL1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COBL like 1;Cordon bleu protein like 1;KIAA977;COBL1;COBLL1 KIAA0977
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cell-cell adherens junction,extracellular exosome,
Tissue Specificity	
Function	similarity:Contains 1 WH2 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night