



NIF3L rabbit pAb

Catalog No	YP-Ab-08974
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	NIF3L1 ALS2CR1 MDS015 My018
Immunogen	Synthesized peptide derived from human NIF3L AA range: 85-135
Specificity	This antibody detects endogenous levels of NIF3L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALS2CR1; Amyotrophic lateral sclerosis 2 (juvenile) chromosome region; candidate 1; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein; CALS-7; MDS15; Ngg1 interacting factor 3-like 1; NIF3 (Ngg1 interacting factor 3; S. pombe homolog)-like 1; NIF3 NGG1 interacting factor 3-like 1 (S. pombe); NIF3 NGG1 interacting factor 3-like 1; NIF3-like protein 1; NIF3L_HUMAN; Nif3l1; OTTHUMP2619; OTTHUMP26191; OTTHUMP26192; NIF3L; NIF3L1 ALS2CR1 MDS015 My018
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Interaction with COPS2 may regulate localization to the nucleus. .
Tissue Specificity	
Function	similarity: Belongs to the UPF0135 (NIF3) family., subunit: Homodimer. Interacts with COPS2 (By similarity). Interacts with THOC7.,
Background	

**matters needing attention**


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night