




ADCL2 rabbit pAb

Catalog No	YP-Ab-08976
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AADACL2
Immunogen	Synthesized peptide derived from human ADCL2 AA range: 74-124
Specificity	This antibody detects endogenous levels of ADCL2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADCL2;ARCL1A;ARMD3;Dance;Developmental arteries and neural crest EGF like protein;Developmental arteries and neural crest EGF-like protein;EVEC;Fbln5;FBLN5_HUMAN;FIBL 5;FIBL-5;Fibulin-5;FLJ959;UP5;Urine p5 protein;AADACL2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	
Function	similarity:Belongs to the 'GDXG' lipolytic enzyme family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night