




GLYL2 rabbit pAb

Catalog No	YP-Ab-08982
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GLYATL2
Immunogen	Synthesized peptide derived from human GLYL2 AA range: 84-134
Specificity	This antibody detects endogenous levels of GLYL2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GLYL2;GLYATL2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Expressed at highest levels in salivary gland and trachea. Also detected in thyroid gland, spinal cord, prostate, lung and fetal brain.
Function	catalytic activity:Acyl-CoA + glycine = CoA + N-acylglycine., function:Mitochondrial acyltransferase which transfers the acyl group to the N-terminus of glycine. Can conjugate a multitude of substrates to form a variety of N-acylglycines.,similarity:Belongs to the glycine N-acyltransferase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night