




SDSL rabbit pAb

Catalog No	YP-Ab-08996
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SDSL
Immunogen	Synthesized peptide derived from human SDSL AA range: 250-300
Specificity	This antibody detects endogenous levels of SDSL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L-serine dehydratase/L-threonine deaminase;L-serine deaminase;L threonine dehydratase;L serine deaminase;EC:4.3.1.19;SDSL
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,mitochondrion,extracellular exosome,
Tissue Specificity	
Function	cofactor:Pyridoxal phosphate.,similarity:Belongs to the serine/threonine dehydratase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night