




EMC3 rabbit pAb

Catalog No	YP-Ab-09006
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	EMC3 TMEM111
Immunogen	Synthesized peptide derived from human EMC3 AA range: 54-104
Specificity	This antibody detects endogenous levels of EMC3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EMC3; POB; EMC3 TMEM111
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	sequence caution: Translated as Trp., similarity: Belongs to the UPF0347 (TMEM111) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night