




CF062 rabbit pAb

Catalog No	YP-Ab-09018
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C6orf62 XTP12 Nbla00237
Immunogen	Synthesized peptide derived from human CF062 AA range: 155-205
Specificity	This antibody detects endogenous levels of CF062 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CF062;C6orf62 XTP12 Nbla00237
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night