




CHUR rabbit pAb

Catalog No	YP-Ab-09019
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CHURC1 C14orf52 CHCH My015
Immunogen	Synthesized peptide derived from human CHUR AA range: 51-101
Specificity	This antibody detects endogenous levels of CHUR at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHUR;CHURC1 C14orf52 CHCH My015
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function:Transcriptional activator that mediates FGF signaling during neural development. Plays a role in the regulation of cell movement (By similarity). Does not bind DNA by itself.,mass spectrometry:The mass difference may be explained by bound zinc ions PubMed:17610897,similarity:Belongs to the Churchill family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night