




PRS37 rabbit pAb

Catalog No	YP-Ab-09020
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PRSS37 TRYX2
Immunogen	Synthesized peptide derived from human PRS37 AA range: 30-80
Specificity	This antibody detects endogenous levels of PRS37 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRS37;PRSS37 TRYX2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Secreted .
Tissue Specificity	Testis-specific (PubMed:27649891). Expressed in spermatids (at protein level) (PubMed:27649891).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night