




CCD53 rabbit pAb

Catalog No	YP-Ab-09026
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CCDC53 AD-016 CGI-116 x0009
Immunogen	Synthesized peptide derived from human CCD53 AA range: 81-131
Specificity	This antibody detects endogenous levels of CCD53 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AD-16;CCD53_HUMAN;CCDC53;CGI-116;Coiled-coil domain-containing protein 53;WASH complex subunit CCDC53;CCD53;CCDC53 AD-016 CGI-116 x0009
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Early endosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night