




TTC29 rabbit pAb

Catalog No	YP-Ab-09028
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TTC29
Immunogen	Synthesized peptide derived from human TTC29 AA range: 213-263
Specificity	This antibody detects endogenous levels of TTC29 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTC29
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, flagellum axoneme .
Tissue Specificity	Expressed in spermatozoa (at protein level).
Function	similarity: Contains 5 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night